

Anti-E-Cadherin/Cd324 Monoclonal Antibody For Flow Cytometry - P12830

Item Number: CM0004652



Introduction

-conjugated rabbit monoclonal antibody against human E-Cadherin/CD324 for flow cytometry. Targets extracellular domain (aa 155-707). Species reactivity: human. Predicted molecular weight ~97 kDa. Suitable for cell adhesion and cancer studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human E-Cadherin/CD324 mAb
Remarks/Alias	UVO; CDHE; ECAD; LCAM; Arc-1; BCDS1; CD324
Species	Human
GeneID (human)	999
GeneID	999
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 155-707 of human E-Cadherin/CD324 (NP_004351.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P12830
Protein Weight	90kDa/97kDa
Shipping	Ice bag