

Anti-E-Selectin/Cd62E Polyclonal Antibody For Wb, Elisa - Q00690

Item Number: CM0022253

Introduction

Rabbit polyclonal antibody to mouse E-Selectin/CD62E. Validated for Western blotting and ELISA. Suitable for studying leukocyte adhesion, inflammation, and endothelial activation.

[Learn More](#)

Parameter	Value
Product Name	E-Selectin/CD62E Rabbit pAb
Remarks/Alias	Elam; CD62E; ELAM-1; LECAM2; E-selectin; E-Selectin/CD62E
Species	Mouse
Gene ID (Human)	6401
Gene ID (Mouse)	20339
Immunogen	Recombinant protein corresponding to a sequence within amino acids 22-557 of mouse E-Selectin/CD62E (NP_035475.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q00690
Protein Weight	66kDa
Shipping	Ice bag