


Anti-E2A/Tcf3 Rabbit Monoclonal Antibody For Wb, Elisa - P15923

Item Number: CM0005971



Introduction

High-quality rabbit monoclonal antibody against human E2A/TCF3, validated for WB and ELISA. Suitable for research on transcription regulation, B-cell differentiation, and cancer. Cross-reacts with human and mouse.

[Learn More](#)

Parameter	Value
Product Name	E2A/TCF3 Rabbit mAb
Remarks/Alias	TCF3; AGM8; E2A; E47; ITF1; TCF-3; VDIR; bHLHb21; transcription factor 3
Species	Human
Gene ID (Human)	6929
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 379-525 of human E2A/TCF3 (NP_003191.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P15923
Protein Weight	67kDa/68kDa
Shipping	Ice bag