


Anti-E2F1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q01094

Item Number: CM0020482



Introduction

Rabbit polyclonal antibody targeting E2F1, a transcription factor regulating cell cycle and apoptosis. Validated for WB, IF/ICC, ELISA, cross-reacts with human, mouse, rat. Ideal for oncology and cell biology research.

[Learn More](#)

Parameter	Value
Product Name	E2F1 Rabbit pAb
Remarks/Alias	RBP3; E2F-1; RBAP1; RBBP3; E2F1
Species	Human
GeneID (human)	1869
GeneID	1869
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-301 of human E2F1 (NP_005216.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01094
Protein Weight	47kDa
Shipping	Ice bag