


Anti-E2F2 Rabbit Monoclonal Antibody For Wb And Elisa - Q14209

Item Number: CM0005707



Introduction

High-quality rabbit monoclonal antibody against transcription factor E2F2, validated for Western blot and ELISA in human and mouse samples. Ideal for cell cycle and cancer research. Ships on ice.

[Learn More](#)

Parameter	Value
Product Name	E2F2 Rabbit mAb
Remarks/Alias	E2F-2; E2F2
Species	Human
GeneID (human)	1870
GeneID	1870
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 303-395 of human E2F2 (Q14209).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q14209
Protein Weight	48kDa
Shipping	Ice bag