


Anti-E2F4 Rabbit Monoclonal Antibody For Wb And Elisa - Q16254

Item Number: CM0003559



Introduction

Explore our anti-E2F4 rabbit monoclonal antibody optimized for WB and ELISA in human and mouse models. Ideal for investigating cell cycle regulation and transcriptional control mechanisms.

[Learn More](#)

Parameter	Value
Product Name	E2F4 Rabbit mAb
Remarks/Alias	E2F-4; E2F4
Species	Human
Gene ID	1874
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 250-350 of human E2F4 (Q16254)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q16254
Protein Weight	44kDa
Shipping	Ice bag