


Anti-E2F8 Polyclonal Antibody For Wb And Elisa - A0Avk6

Item Number: CM0015090



Introduction

Rabbit polyclonal antibody against human E2F8, suitable for Western blot and ELISA applications. Cross-reacts with mouse. Validated for research on cell cycle, transcription regulation, and angiogenesis.

[Learn More](#)

Parameter	Value
Product Name	E2F8 Rabbit pAb
Remarks/Alias	E2F-8; E2F8
Species	Human
GeneID (Human)	79733
GeneID	79733
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 788-867 of human E2F8 (NP_078956.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	A0AVK6
Protein Weight	94kDa
Shipping	Ice bag