


Anti-Eaat1/Slc1A3 Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P43003

Item Number: CM0005614



Introduction

Rabbit monoclonal antibody to human EAAT1/SLC1A3, validated for WB, IF-P, IHC-P, ELISA. Recognizes mouse and rat. Suitable for neuroscience studies on glutamate transport and excitatory neurotransmission.

[Learn More](#)

Parameter	Value
Product Name	EAAT1/SLC1A3 Rabbit mAb
Remarks/Alias	EA6; EAAT1; GLAST; GLAST1; EAAT1/SLC1A3
Species	Human
GeneID (Human)	6507
GeneID	6507
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 150-250 of human EAAT1/SLC1A3 (P43003)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P43003
Protein Weight	60kDa
Shipping	Ice bag