


Anti-Eaat2 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P43004

Item Number: CM0009974



Introduction

High-affinity rabbit monoclonal antibody targeting human EAAT2 (SLC1A2) for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for glutamate transporter studies in synaptic function and neurology research.

[Learn More](#)

Parameter	Value
Product Name	EAAT2 Rabbit mAb
Remarks/Alias	GLT1; HBGT; DEE41; EAAT2; GLT-1; EIEE41; EAAT2/SLC1A2
Species	Human
Gene ID	6506
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43004
Protein Weight	62kDa
Shipping	Ice bag