


Anti-Eaat2/Slc1A2 Polyclonal Antibody For Wb, If-P, Elisa - P43004

Item Number: CM0014539



Introduction

Rabbit polyclonal antibody against EAAT2/SLC1A2 (GLT-1), validated for WB, IF-P, ELISA; reacts with human, mouse, rat. Targets the glutamate transporter for synaptic biology and neurological disease research.

[Learn More](#)

Parameter	Value
Product Name	EAAT2/SLC1A2 Rabbit pAb
Remarks/Alias	GLT1; HBGT; DEE41; EAAT2; GLT-1; EIEE41; EAAT2/SLC1A2
Species	Human
GeneID (Human)	6506
GeneID	6506
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 474-574 of human EAAT2/SLC1A2 (NP_004162.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P43004
Protein Weight	62kDa
Shipping	Ice bag