


Anti-Eaf1 Polyclonal Antibody For Wb, Elisa - Q96Jc9

Item Number: CM0022219



Introduction

Rabbit polyclonal antibody targeting human EAF1 (ELL-associated factor 1), validated for Western blot and ELISA. Cross-reacts with human and rat. EAF1 is a transcriptional transactivator of ELL/ELL2 elongation activities, localized in nuclear speckles and Cajal bodies.

[Learn More](#)

Parameter	Value
Product Name	EAF1 Rabbit pAb
Remarks/Alias	EAF1
Species	Human
GeneID (Human)	85403
GeneID	85403
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-120 of human EAF1 (NP_149074.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q96JC9
Protein Weight	29kDa
Shipping	Ice bag