


Anti-Eb3/Mapre3 Polyclonal Antibody For Wb, Elisa - Q9Upy8

Item Number: CM0029570



Introduction

Rabbit polyclonal antibody targeting human EB3/MAPRE3, validated for Western blotting and ELISA. Cross-reacts with mouse. Designed for microtubule dynamics, spindle function, and cell division studies.

[Learn More](#)

Parameter	Value
Product Name	EB3/MAPRE3 Rabbit pAb
Remarks/Alias	EB3; RP3; EBF3; EBF3-S; EB3/MAPRE3
Species	Human
GeneID (Human)	22924
GeneID	22924
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-281 of human EB3/MAPRE3 (NP_036458.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9UPY8
Protein Weight	32kDa
Shipping	Ice bag