

Anti-Ebag9 Rabbit Polyclonal Antibody For Wb, Elisa - 000559

Item Number: CM0018835

Introduction

Rabbit polyclonal antibody against human EBAG9 (RCAS1), validated for WB and ELISA. Suitable for apoptosis and cancer research, targeting a 24 kDa Golgi membrane protein.

[Learn More](#)

Parameter	Value
Product Name	EBAG9 Rabbit pAb
Remarks/Alias	EB9; PDAF; EBAG9
Species	Human
GeneID (Human)	9166
GeneID	9166
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-213 of human EBAG9 (NP_004206.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O00559
Protein Weight	24kDa
Shipping	Ice bag