

Anti-Ebf3 Polyclonal Antibody For Wb, Elisa - Q9H4W6

Item Number: CM0019070



Introduction

Rabbit polyclonal antibody targeting human EBF3 (COE3). Validated for Western blot and ELISA, with cross-reactivity to rat. Ideal for transcription regulation and neurodevelopmental research.

[Learn More](#)

Parameter	Value
Product Name	EBF3 Rabbit pAb
Remarks/Alias	COE3, OE-2, EBF-3, HADD5, O/E-2, EBF3
Species	Human
GeneID (human)	253738
GeneID	253738
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 400-500 of human EBF3 (NP_001005463.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q9H4W6
Protein Weight	65kDa
Shipping	Ice bag