

Anti-Ebi3 Rabbit Monoclonal Antibody For Wb, Elisa - Q14213

Item Number: CM0002141



Introduction

Rabbit monoclonal antibody against human EBI3 (IL-27 subunit beta), cross-reactive with human, mouse, rat. Validated for WB and ELISA. Target is a secreted cytokine involved in immune regulation.

[Learn More](#)

Parameter	Value
Product Name	EBI3 Rabbit mAb
Remarks/Alias	IL27B; IL35B; IL-27B; EBI3
Species	Human
GeneID (human)	10148
GeneID	10148
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 130-229 of human EBI3 (Q14213)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14213
Protein Weight	25kDa
Shipping	Ice bag