


Anti-Ebp Polyclonal Antibody For Wb And Elisa - Q15125

Item Number: CM0016211



Introduction

High-quality rabbit polyclonal antibody against human EBP (3-beta-hydroxysteroid-Delta(8),Delta(7)-isomerase), validated for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for cholesterol metabolism research.

[Learn More](#)

Parameter	Specification
Product Name	EBP Rabbit pAb
Remarks/Alias	CPX; CHO2; CPXD; MEND; CDPX2; EBP
Species	Human
Gene ID	10682
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 50-150 of human EBP (NP_006570.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q15125
Protein Weight	26kDa
Shipping	Ice bag