

Anti-Ecd Rabbit Polyclonal Antibody For Wb And Elisa - O95905

Item Number: CM0022236

Introduction

High-quality rabbit polyclonal antibody targeting ECD (ecdysoneless homolog) for WB and ELISA. Cross-reacts with human, mouse, and rat. Ideal for studying p53 stability, transcription, and cell cycle progression.

[Learn More](#)

Parameter	Value
Product Name	ECD Rabbit pAb
Remarks/Alias	GCR2; SGT1; HSGT1; ECD
Species	Human
GeneID (Human)	11319
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 435-644 of human ECD (NP_009196.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95905
Protein Weight	73kDa
Shipping	Ice bag