


Anti-Ech1 Rabbit Monoclonal Antibody (Kd Validated) For Wb, Ihc-P, If/Icc, Elisa - Q13011

Item Number: CM0011404



Introduction

KD-validated rabbit monoclonal anti-ECH1 antibody (HPXEL) for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Targets mitochondrial delta(3,5)-delta(2,4)-dienoyl-CoA isomerase (SWISS Q13011).

[Learn More](#)

Parameter	Specification
Product Name	[KD Validated] ECH1 Rabbit mAb
Remarks/Alias	HPXEL
Species	Human
GeneID (Human)	1891
GeneID	1891
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-322 of human ECH1 (NP_001389.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13011
Protein Weight	36kDa
Shipping	Ice bag