

Anti-Echdc2 Rabbit Polyclonal Antibody For Wb And Elisa - Q86Yb7

Item Number: CM0016343

Introduction

ECHDC2 Rabbit Polyclonal Antibody validated for Western blotting (WB) and ELISA. Recognizes human ECHDC2 and cross-reacts with mouse and rat. Suitable for mitochondrial fatty acid and lipid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	ECHDC2 Rabbit pAb
Remarks/Alias	ECHDC2
Species	Human
GeneID (Human)	55268
GeneID	55268
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-292 of human ECHDC2 (NP_001185890.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q86YB7
Protein Weight	31kDa
Shipping	Ice bag