


Anti-Eci1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P42126

Item Number: CM0015165



Introduction

Rabbit polyclonal antibody targeting ECI1 (Enoyl-CoA delta isomerase 1, mitochondrial), validated for WB, IHC-P, IF/ICC, and ELISA; cross-reacts with human and mouse; 33 kDa; ideal for fatty acid metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	ECI1 Rabbit pAb
Remarks/Alias	DCI; ECI1
Species	Human
GeneID (human)	1632
GeneID	1632
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P42126
Protein Weight	33kDa
Shipping	Ice bag