


## Anti-Eci2 Rabbit Polyclonal Antibody For Wb And Elisa - 075521

Item Number: CM0013443



### Introduction

High-quality rabbit polyclonal antibody against human ECI2 (Enoyl-CoA delta isomerase 2). Validated for WB, ELISA with cross-reactivity to mouse and rat. Ideal for research on fatty acid metabolism and peroxisomal disorders.

[Learn More](#)

| Parameter        | Value                                                                                                           |
|------------------|-----------------------------------------------------------------------------------------------------------------|
| Product Name     | ECI2 Rabbit pAb                                                                                                 |
| Remarks/Alias    | DRS1; PECI; ACBD2; DRS-1; HCA88; dj1013A10.3; ECI2                                                              |
| Species          | Human                                                                                                           |
| GeneID (Human)   | 10455                                                                                                           |
| GeneID           | 10455                                                                                                           |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-364 of human ECI2 (NP_006108.2) |
| Source           | Rabbit                                                                                                          |
| Category         | Polyclonal Antibodies                                                                                           |
| Application      | WB, ELISA                                                                                                       |
| Cross Reactivity | Human, Mouse, Rat                                                                                               |
| SWISS            | 075521                                                                                                          |
| Protein Weight   | 44kDa                                                                                                           |
| Shipping         | Ice bag                                                                                                         |