

Anti-Ecm1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q16610

Item Number: CM0022996

Introduction

Rabbit polyclonal antibody against human ECM1, validated for WB, IF/ICC, and ELISA. It cross-reacts with mouse and targets a 61kDa secreted glycoprotein involved in bone formation and angiogenesis, suitable for extracellular matrix and cancer research.

[Learn More](#)

Parameter	Value
Product Name	ECM1 Rabbit pAb
Remarks/Alias	URBWD; ECM1
Species	Human
GeneID (Human)	1893
GeneID	1893
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 280-405 of human ECM1 (NP_004416.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q16610
Protein Weight	61kDa
Shipping	Ice bag