

Anti-Ecsit Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Bq95

Item Number: CM0028982



Introduction

Rabbit polyclonal antibody targeting ECSIT (SITPEC), validated for WB, IF/ICC, and ELISA in human, mouse, and rat. ECSIT is a key adaptor in innate immunity, NF-κB signaling, and mitochondrial complex I assembly.

[Learn More](#)

Parameter	Value
Product Name	ECSIT Rabbit pAb
Remarks/Alias	SITPEC; ECSIT
Species	Human
Gene ID (Human)	51295
Gene ID	51295
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-296 of human ECSIT (NP_001135936.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BQ95
Protein Weight	49kDa
Shipping	Ice bag