

Anti-Edar Rabbit Polyclonal Antibody For Wb, Elisa - Q9Une0

Item Number: CM0028085



Introduction

Rabbit polyclonal antibody targeting human EDAR, the ectodysplasin A receptor. Suitable for Western blot and ELISA applications. EDAR mediates NF-kB and JNK signaling and is implicated in ectodermal dysplasia.

[Learn More](#)

Parameter	Value
Product Name	EDAR Rabbit pAb
Remarks/Alias	DL; ED3; ED5; ED1R; EDA3; HRM1; EDA1R; ECTD10A; ECTD10B; EDA-A1R; EDAR
Species	Human
GeneID (human)	10913
GeneID	10913
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-187 of human EDAR (NP_071731.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UNE0
Protein Weight	49kDa
Shipping	Ice bag