

Anti-Edaradd Rabbit Polyclonal Antibody For Wb, Elisa - Q8Wwz3

Item Number: CM0022684

Introduction

Rabbit polyclonal EDARADD antibody for WB and ELISA applications. Validated in human, cross-reacts with mouse and rat. Ideal for ectodermal dysplasia and NF-kB signaling research.

[Learn More](#)

Parameter	Value
Product Name	EDARADD Rabbit pAb
Remarks/Alias	CR; ED3; EDA3; ECTD11A; ECTD11B; EDARADD
Species	Human
GeneID (Human)	128178
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-215 of human EDARADD (NP_665860.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8WWZ3
Protein Weight	25kDa
Shipping	Ice bag