

Anti-Edc3 Rabbit Polyclonal Antibody For Wb - Q96F86

Item Number: CM0013059



Introduction

EDC3 rabbit polyclonal antibody validated for WB and ELISA. Detects human, mouse, rat EDC3 (56kDa). Recombinant immunogen (aa159-508) ensures specificity. Ideal for mRNA decapping, P-body, and spermiogenesis research.

[Learn More](#)

Parameter	Value
Product Name	EDC3 Rabbit pAb
Remarks/Alias	YJDC; LSM16; MRT50; YJEFN2; hYjeF_N2-15q23; EDC3
Species	Human
Gene ID (Human)	80153
Gene ID	80153
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 159-508 of human EDC3 (NP_079359.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96F86
Protein Weight	56kDa
Shipping	Ice bag