

Anti-Edem2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Bv94

Item Number: CM0022970

Introduction

Rabbit polyclonal antibody against human EDEM2, an ERAD pathway component. Validated for WB, IF/ICC, and ELISA; cross-reactive with mouse. For studying ER-associated degradation and glycoprotein quality control.

[Learn More](#)

Parameter	Value
Product Name	EDEM2 Rabbit pAb
Remarks/Alias	C20orf31; C20orf49; bA4204.1; EDEM2
Species	Human
GenelD (human)	55741
GenelD	55741
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 399-578 of human EDEM2 (NP_060687.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BV94
Protein Weight	65kDa
Shipping	Ice bag