


Anti-Edf1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O60869

Item Number: CM0019948



Introduction

Rabbit polyclonal antibody targeting human EDF1, validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for endothelial differentiation and transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] EDF1 Rabbit pAb
Remarks/Alias	MBF1; EDF-1; CFAP280; F1
Species	Human
GeneID (Human)	8721
GeneID	8721
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-148 of human EDF1 (NP_003783.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60869
Protein Weight	16kDa
Shipping	Ice bag