

## Anti-Edg1/Hexim1 Monoclonal Antibody For Wb And Elisa - P21453

Item Number: CM0008606



### Introduction

Rabbit monoclonal antibody targeting EDG1/HEXIM1 (S1PR1), suitable for WB and ELISA with cross-reactivity to human, mouse, and rat. Ideal for GPCR signaling and angiogenesis research.

[Learn More](#)

| Parameter        | Value                                                                                                                      |
|------------------|----------------------------------------------------------------------------------------------------------------------------|
| Product Name     | EDG1/HEXIM1 Rabbit mAb                                                                                                     |
| Remarks/Alias    | EDG1; S1P1; CD363; ECGF1; EDG-1; CHEDG1; D153362; EDG1/HEXIM1                                                              |
| Species          | Human                                                                                                                      |
| Gene ID (Human)  | 1901                                                                                                                       |
| Gene ID          | 1901                                                                                                                       |
| Immunogen        | Synthetic Peptide  A synthetic peptide corresponding to a sequence within amino acids 1-100 of human EDG1/HEXIM1 (P21453). |
| Source           | Rabbit                                                                                                                     |
| Category         | Monoclonal Antibodies                                                                                                      |
| Application      | WB, ELISA                                                                                                                  |
| Cross Reactivity | Human, Mouse, Rat                                                                                                          |
| SWISS            | P21453                                                                                                                     |
| Protein Weight   | 43kDa                                                                                                                      |
| Shipping         | Ice bag                                                                                                                    |