

Anti-Edg1/S1Pr1 Polyclonal Antibody For Wb And Elisa - P21453

Item Number: CM0018654



Introduction

High-quality rabbit polyclonal antibody targeting human S1PR1 (EDG1, S1P1). Validated for WB and ELISA; cross-reacts with mouse. Ideal for studying GPCR signaling in cell migration, angiogenesis, and immune response.

[Learn More](#)

Parameter	Value
Product Name	EDG1/HEXIM1 Rabbit pAb
Remarks/Alias	EDG1; S1P1; CD363; ECGF1; EDG-1; CHEDG1; D1S3362; EDG1/HEXIM1
Species	Human
GeneID (Human)	1901
GeneID	1901
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 303-382 of human EDG1/HEXIM1 (NP_001391.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P21453
Protein Weight	43kDa
Shipping	Ice bag