

Anti-Edg2/Lpa-1 (Lysophosphatidic Acid Receptor 1) Rabbit Monoclonal Antibody For Wb, Elisa - Q92633

Item Number: CM0002879



Introduction

High-quality rabbit monoclonal anti-EDG2/LPA-1 (LPAR1) antibody for Western blot (WB) and ELISA applications. Validated in human samples, reacts with human EDG2/LPA-1. Suitable for cell signaling, migration, and receptor research.

[Learn More](#)

Parameter	Value
Product Name	EDG2/LPA-1 Rabbit mAb
Remarks/Alias	EDG2; LPA1; VZG1; edg-2; vzg-1; Gpcr26; Mrec1.3; rec.1.3; EDG2/LPA-1
Species	Human
Gene ID (Human)	1902
Gene ID	1902
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human EDG2/LPA-1 (Q92633).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q92633
Protein Weight	41kDa
Shipping	Ice bag