


Anti-Edn1 Polyclonal Antibody (Ko Validated) For Wb, If/Icc, Elisa - P05305

Item Number: CM0021365



Introduction

EDN1 Rabbit Polyclonal Antibody, KO Validated, for Western Blot, IF/ICC, ELISA. Reacts with Human, Mouse, Rat. Targets Endothelin-1, a vasoconstrictor peptide. Ideal for IVD research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] EDN1 Rabbit pAb
Remarks/Alias	ET1; QME; PPET1; ARCND3; HDLCQ7; N1
Species	Human
GeneID (Human)	1906
GeneID	1906
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-212 of human EDN1 (NP_001946.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05305
Protein Weight	24kDa
Shipping	Ice bag