

Anti-Ednrb Polyclonal Antibody For Wb And Elisa - P24530

Item Number: CM0027117



Introduction

Rabbit polyclonal antibody targeting human EDNRB, validated for WB and ELISA. Cross-reacts with human and mouse. Ideal for endothelin receptor signaling, vascular research, and developmental studies.

[Learn More](#)

Parameter	Value
Product Name	EDNRB Rabbit pAb
Remarks/Alias	ETB; ET-B; ETB1; ETBR; ETRB; HSCR; WS4A; ABCDS; ET-BR; HSCR2; EDNRB
Species	Human
GeneID	1910
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human EDNRB (NP_001116131.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P24530
Protein Weight	50kDa
Shipping	Ice bag