


Anti-Eea1 Rabbit Monoclonal Antibody For Wb, Elisa - Q15075

Item Number: CM0011510



Introduction

Rabbit monoclonal antibody against human EEA1, suitable for WB and ELISA with cross-reactivity to mouse and rat. The ~162 kDa early endosome marker is essential for endosomal trafficking studies.

[Learn More](#)

Parameter	Value
Product Name	EEA1 Rabbit mAb
Remarks/Alias	MST105; ZFYVE2; MSTP105; EEA1
Species	Human
Gene ID	8411
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human EEA1 (Q15075)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15075
Protein Weight	162kDa
Shipping	Ice bag