


Anti-Eea1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q15075

Item Number: CM0016168



Introduction

Rabbit polyclonal antibody targeting human EEA1, validated in human, mouse, and rat samples. Suitable for Western blot, immunohistochemistry, and ELISA applications in endosomal trafficking research.

[Learn More](#)

Parameter	Value
Product Name	EEA1 Rabbit pAb
Remarks/Alias	MST105; ZFYVE2; MSTP105; EEA1
Species	Human
Gene ID (Human)	8411
Gene ID	8411
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1182-1411 of human EEA1 (NP_003557.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15075
Protein Weight	162kDa
Shipping	Ice bag