

Anti-Eed Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O75530

Item Number: CM0014414



Introduction

Rabbit polyclonal antibody targeting human EED protein. Validated for Western blot, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for chromatin regulation and cancer research.

[Learn More](#)

Parameter	Value
Product Name	EED Rabbit pAb
Remarks/Alias	HEED; COGIS; WAIT1; EED
Species	Human
GeneID (Human)	8726
Immunogen	Synthetic peptide within amino acids 1-100 of human EED (NP_003788.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75530
Protein Weight	50kDa
Shipping	Ice bag