


Anti-Eef1A1 Polyclonal Antibody For Wb, If/Icc, Elisa - P68104

Item Number: CM0018645



Introduction

eEF1A1 Rabbit pAb validated for WB, IF/ICC, ELISA, reacting with human, mouse, rat. This polyclonal antibody targets translation elongation factor eEF1A1 (P68104), suitable for research on protein synthesis and host-virus interactions.

[Learn More](#)

Parameter	Value
Product Name	eEF1A1 Rabbit pAb
Remarks/Alias	CCS3; EF1A; PTI1; CCS-3; EE1A1; EEF-1; EEF1A; EF-Tu; EF1A1; LENG7; eEF1A-1; GRAF-1EF; EF1alpha1; eEF1A1
Species	Human
GeneID (human)	1915
GeneID	1915
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 123-462 of human eEF1A1 (NP_001393.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P68104
Protein Weight	50kDa
Shipping	Ice bag