


Anti-Eef1A1 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P68104

Item Number: CM0007211



Introduction

Rabbit monoclonal antibody against human eEF1A1, validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for protein synthesis and translation research.

[Learn More](#)

Parameter	Value
Product Name	eEF1A1 Rabbit mAb
Remarks/Alias	CCS3; EF1A; PTI1; CCS-3; EE1A1; EEF-1; EEF1A; EF-Tu; EF1A1; LENG7; eEF1A-1; GRAF-1EF; EF1alpha1; eEF1A1
Species	Human
Gene ID	1915
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-400 of human eEF1A1. (NP_001393.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P68104
Protein Weight	50kDa
Shipping	Ice bag