


Anti-Eef1A1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P68104

Item Number: CM0020998



Introduction

Rabbit polyclonal antibody targeting human eEF1A1 (Elongation factor 1-alpha 1). Suitable for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Useful for translation elongation and protein synthesis research.

[Learn More](#)

Parameter	Specification
Product Name	eEF1A1 Rabbit pAb
Remarks/Alias	CCS3; EF1A; PTI1; CCS-3; EE1A1; EEF-1; EEF1A; EF-Tu; EF1A1; LENG7; eEF1A-1; GRAF-1EF; EF1alpha1; eEF1A1
Species	Human
Gene ID	1915
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 301-400 of human eEF1A1 (NP_001393.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P68104
Protein Weight	50kDa
Shipping	Ice bag