


Anti-Eef1D Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P29692

Item Number: CM0018764



Introduction

Rabbit polyclonal antibody targeting elongation factor 1-delta (EEF1D), validated for WB, IHC-P, IF/ICC, IP, ELISA. Reacts with human, mouse, rat. Relevant for protein synthesis, transcription regulation, and heat shock response research.

[Learn More](#)

Parameter	Specification
Product Name	EEF1D Rabbit pAb
Remarks/Alias	EF1D; EF-1D; FP1047; EEF1D
Species	Human
GeneID (Human)	1936
GeneID	1936
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human EEF1D (NP_001123529.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29692
Protein Weight	31kDa
Shipping	Ice bag