


Anti-Eef1E1 Polyclonal Antibody For Wb, If/Icc, Elisa - O43324

Item Number: CM0013093



Introduction

High-specificity rabbit polyclonal antibody against human EEF1E1 (AIMP3/p18), validated for WB, IF/ICC, and ELISA. Detects human, mouse, and rat targets. Ideal for DNA damage response and cancer research.

[Learn More](#)

Parameter	Value
Product Name	EEF1E1 Rabbit pAb
Remarks/Alias	P18; AIMP3; EEF1E1
Species	Human
GeneID (Human)	9521
GeneID	9521
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human EEF1E1 (NP_004271.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43324
Protein Weight	20kDa
Shipping	Ice bag