


Anti-Eef1G Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P26641

Item Number: CM0006053



Introduction

Rabbit monoclonal antibody against human elongation factor 1-gamma (eEF1G). Validated for WB, IHC-P, IF/ICC, and ELISA; cross-reacts with mouse and rat. For protein biosynthesis and cancer research.

[Learn More](#)

Parameter	Value
Product Name	eEF1G Rabbit mAb
Remarks/Alias	EF1G; GIG35; eEF1G
Species	Human
Gene ID (Human)	1937
Gene ID	1937
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-123 of human eEF1G (P26641).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26641
Protein Weight	50kDa
Shipping	Ice bag