


Anti-Eef1G Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P26641

Item Number: CM0030935



Introduction

Rabbit polyclonal antibody targeting eEF1G (EF1G, GIG35) for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Suitable for research on protein biosynthesis and elongation factor complex.

[Learn More](#)

Parameter	Value
Product Name	eEF1G Rabbit pAb
Remarks/Alias	EF1G; GIG35; eEF1G
Species	Human
GeneID (Human)	1937
GeneID	1937
Immunogen	Recombinant protein (sequence: amino acids 168-437 of human eEF1G, NP_001395.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26641
Protein Weight	50kDa
Shipping	Ice bag