


Anti-Eef2 Rabbit Polyclonal Antibody For Wb And Elisa - P13639

Item Number: CM0022353



Introduction

Rabbit polyclonal antibody targeting human elongation factor 2 (eEF2) for Western blot and ELISA. Reacts with human eEF2; immunogen recombinant protein (aa 1-290). Ideal for translation elongation studies and SCA26 research. Predicted molecular weight 95 kDa.

[Learn More](#)

Parameter	Value
Product Name	eEF2 Rabbit pAb
Remarks/Alias	EF2; EF-2; EEF-2; SCA26; eEF2
Species	Human
Gene ID (Human)	1938
Gene ID	1938
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human eEF2 (NP_001952.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P13639
Protein Weight	95kDa
Shipping	Ice bag