


Anti-Eef2K Polyclonal Antibody For Wb, If/Icc, Elisa - 000418

Item Number: CM0025344



Introduction

Validated Anti-EEF2K rabbit polyclonal antibody for WB, IF/ICC, ELISA. Recognizes human and mouse EEF2K; 82kDa. Ideal for studying translation regulation and kinase signaling.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] EEF2K Rabbit pAb
Remarks/Alias	eEF-2K; CaMKIII; HSU93850; 2K
Species	Human
GeneID (Human)	29904
GeneID	29904
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 406-725 of human EEF2K (NP_037434.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O00418
Protein Weight	82kDa
Shipping	Ice bag