

Anti-Efcab3 Rabbit Polyclonal Antibody For Wb, Elisa - Q8N7B9

Item Number: CM0014682

Introduction

EFCAB3 Rabbit Polyclonal Antibody (pAb) raised against recombinant human EFCAB3. Validated for WB and ELISA. Reacts with human, mouse, and rat samples. Suitable for calcium-binding domain studies.

[Learn More](#)

Parameter	Value
Product Name	EFCAB3 Rabbit pAb
Remarks/Alias	EFCAB3
Species	Human
GeneID (Human)	146779
GeneID	146779
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 249-438 of human EFCAB3 (NP_775774.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8N7B9
Protein Weight	50kDa
Shipping	Ice bag