

Anti-Efhc1 Rabbit Polyclonal Antibody For Wb, Elisa - Q5JvL4

Item Number: CM0029073



Introduction

Anti-EFHC1 Rabbit Polyclonal Antibody detects EFHC1 in human, mouse, rat samples by Western blot and ELISA. Useful for cilia, neuronal migration, and epilepsy research.

[Learn More](#)

Parameter	Value
Product Name	EFHC1 Rabbit pAb
Remarks/Alias	EJM1, POC9, RIB72, dj304B14.2, EFHC1
Species	Human
GeneID (Human)	114327
GeneID	114327
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 391-640 of human EFHC1 (NP_060570.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5JVL4
Protein Weight	74kDa
Shipping	Ice bag