

Anti-Efhd2 Polyclonal Antibody For Wb, Elisa - Q96C19

Item Number: CM0015380



Introduction

Rabbit polyclonal antibody targeting EFHD2 (Swiprosin-1) for Western blot and ELISA. Reacts with human and mouse samples; useful for studying B-cell apoptosis and NF-kappa-B regulation.

[Learn More](#)

Parameter	Value
Product Name	EFHD2 Rabbit pAb
Remarks/Alias	SWS1; EFHD2
Species	Human
Gene ID (Human)	79180
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-78 of human EFHD2 (NP_077305.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q96C19
Protein Weight	27kDa
Shipping	Ice bag