


Anti-Efna1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P20827

Item Number: CM0026014



Introduction

Rabbit polyclonal antibody targeting human EFNA1 (ephrin-A1). Validated for WB, IHC-P, and ELISA; cross-reacts with mouse EFNA1. Suitable for studying angiogenesis, tumor biology, and Eph receptor signaling.

[Learn More](#)

Parameter	Value
Product Name	EFNA1 Rabbit pAb
Remarks/Alias	B61; EFL1; GMAN; ECKLG; EPLG1; LERK1; LERK-1; TNFAIP4; EFNA1
Species	Human
GeneID (Human)	1942
GeneID	1942
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 19-182 of human EFNA1 (NP_004419.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P20827
Protein Weight	24kDa
Shipping	Ice bag