

Anti-Efna3 Rabbit Polyclonal Antibody For Wb, Elisa - P52797

Item Number: CM0028518

Introduction

Anti-EFNA3 rabbit polyclonal antibody validated for WB and ELISA with mouse cross-reactivity. Targets human ephrin-A3, a GPI-anchored Eph receptor ligand crucial in neuronal and vascular development. Ideal for IVD development.

[Learn More](#)

Parameter	Value
Product Name	EFNA3 Rabbit pAb
Remarks/Alias	EFL2; EPLG3; LERK3; Ehk1-L; EFNA3
Species	Human
GeneID (Human)	1944
GeneID	1944
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 23-214 of human EFNA3 (NP_004943.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P52797
Protein Weight	26kDa
Shipping	Ice bag